

Optimizing Spectral Differentiation in Super-Resolution Imaging via Chernoff Information

Matthew Hogan, MS (Presenting)
Sanduni Fernando, PhD
Jordan Gerton, PhD
Erik Jorgensen, PhD

March 19, 2026




Global Physics
Summit



The Diffraction Limit

Green Fluorescent Protein (GFP) is small



Even smaller than the wavelength of light it emits (~500 nm)


The Diffraction Limit

Interference results in a point-spread function (PSF)

For a perfectly circular aperture, this is known as the Airy disk

The Diffraction Limit


Structures smaller than the critical distance cannot
be resolved by traditional light microscopes



Ernst Abbe (1873): $d \geq \frac{\lambda}{2NA} \approx 200 \text{ nm}$


Single Molecule Localization Microscopy (SMLM)

Measure




(Disclaimer: PSF is simulated)

Fit



Reconstruct




STORM, PALM, dna-PAINT, etc


Spatial resolution of the reconstructed image is, at best, the *Cramer-Rao Bound* (~ 10 nm)

What about color information?

Simultaneous Multicolor SMLM?




PSF Engineering with a Glass Phase Plate




It Really Works!

Good distinguishability for 2 colors



Worse distinguishability for 3 colors




What's the best we could do?

Chernoff Information


Suppose we perform hypothesis tests to classify

$$\nu =$$




as either

$$\nu_1 =$$



or

$$\nu_2 =$$



The Chernoff Information between an observation ν and a known class ν_i (under some assumptions) is

$$C_i(\beta) = \nu_i \frac{(\beta - 1)(\log(\frac{\beta - 1}{\log \beta}) - 1) + \log \beta}{\log \beta}; \quad \beta = \frac{\nu}{\nu_i}$$


If $C_1 < C_2 \implies \nu \rightarrow \nu_1$, else $\nu \rightarrow \nu_2$

Chernoff Information tightly bounds Bayes' probability of error when classifying via hypothesis testing (similar to CRB/Fisher Information)


Tightest bound = best decision
(least probability of error)

Symmetric X-PP is Optimal for 2 Colors



Original X-PP





(Optimal for PSFs of wavelengths $\lambda = 440 \text{ nm}$ and $\lambda = 880 \text{ nm}$)



Chernoff-Information-Maximized X-PP



Symmetric Hex-PP is Optimal for 3 Colors



Diffraction-Limited Fluorescent Microspheres

The **0.1 μm TetraSpeck™ microspheres** are stained throughout with four different fluorescent dyes, yielding beads that each display four well-separated excitation/emission peaks - 360/430 nm (blue), 505/515 nm (green), 560/580 nm (orange) and 660/680 (dark red).




$$\lambda_{exc} = 488 \text{ nm}$$

$$\lambda_{exc} = 568 \text{ nm}$$

$$\lambda_{exc} = 647 \text{ nm}$$

Improved Spectral Classification...

... at the Cost of Spatial Resolution



Thank You!